

3M Personal Safety Division



Design of Filtering Facepiece Respirators

NIOSH Healthcare Stakeholder Meeting

Lauri Alvarez – 3M Company

June 18, 2013



Design of Filtering Facepiece Respirators

- 3M Involvement in Project BREATHE
- 3M Technologies and Testing Capabilities related to development of health care respirators



3M and Project BREATHE

- Kick-off meeting to work with VA and NIOSH March 2012
- 3M Generated Concepts – July 2012

23 Ideas Related To:

Noseclip

Headbands



Respirator Shape
(Style)

Adhesive seal

Nosefoam

Resulting in 12 Concepts Proposed to VA



Concept Evaluation

- 12 Concepts to VA for Idea Screening Compared to 3M 1860

| Question | Much Worse | Worse | The Same | Better | Much Better |
|--|------------|-------|----------|--------|-------------|
| Overall, my opinion is the Prototype Respirator compared to 1860 is: | | | | | |
| Overall, the comfort of the Prototype Respirator compared to 1860 is: | | | | | |
| Ease of donning (putting on) of Prototype compared to 1860 is: | | | | | |
| Ease of doffing (removing) of Prototype compared to 1860 is: | | | | | |
| Comfort of Prototype respirator (facial heat) compared to 1860 is: | | | | | |
| Comfort in the nose area of Prototype respirator compared to 1860 is: | | | | | |
| Comfort on face (scratching, itching) of prototype compared to 1860: | | | | | |
| The ease of breathing through prototype respirator compared to 1860: | | | | | |
| How well does this prototype hold its shape compared to 1860 | | | | | |
| Weight of the prototype on the face compared to 1860 | | | | | |
| How do headbands of prototype respirator compare to 1860? Explain in comment area below in detail. | | | | | |
| Field of view of prototype respirator compared to 1860 is: | | | | | |



Concept to Prototype


- Idea Screening by VA Participants (Office Setting)

12 Concepts

Comfort on
Nose

Shape Retention

Opinion of Headbands



Ease of Breathing

Facial Heat

Overall Opinion

Resulting in 2 Prototypes –
1 Cup Style, 1 Flat Fold

Next Steps -

- 3M will provide 2 prototypes
- VA will conduct Simulated Workplace Study
 - Evaluate Fit and Comfort
 - Tasks to Represent Patient Care Activities



3M Technologies and Testing Capabilities

- Respirator Fit and Comfort
 - Fit Testing
 - Thermal Imaging
- Filtration Technology
- Biological aerosol interaction with respirators
- Fluid Resistance



Respirator Fit

- Consensus 4: Respirator Fit
- Consensus 9: Gauging Fit
- 3M has internal fit testing capabilities
- We have an ongoing technology program to develop test methods to evaluate respirator fit and comfort



3M Personal Safety Division

3M Fit Testing Laboratory

- Multiple fit test chambers



Simulated Workplace Study in Health Care Setting

- Collaborative project with Professor Lisa Brosseau (University of Minnesota, School of Public Health)
- Objective:
 - To measure simulated workplace fit factors for filtering facepiece respirators during typical patient care tasks (e.g. patient assessment, IV care, wound care) conducted by healthcare workers in a simulated healthcare environment using facilities at the Interprofessional Education and Research Center (IERC), University of Minnesota



3M Personal Safety Division


Subjects: 8 registered nurses

Respirators: 1860 and 1860S

Continuous fit monitoring to calculate SWFF



Scenario Geometric Mean SWFF and FF by Subject



- Simulated Workplace Fit Factor (SWFF) highly dependent on the individual subject.
- Overall, SWFF is higher than Fit Test calculated Fit Factor (FF).

Jessica Hauge, Marc Roe, Lisa M. Brosseau & Craig Colton (2012): Real-Time Fit of a Respirator during Simulated Health Care Tasks, Journal of Occupational and Environmental Hygiene, 9:10, 563-571



Respirator Comfort



- Consensus 19: Facial Heat
- 3M developed real-time thermal imaging methods to identify effect of respirator design on heat build up



Improving Respirator Comfort – Valve Performance

Test method developed to identify effect of valve on heat build up in respirator


- Red heads are supplied with humidified, body temperature breathing air
- External temperature measured using thermal imaging camera
- Temperature inside the respirator and flow resistance measured using sensors embedded in the red head



3M Personal Safety Division


Filtration Technology

- Electrostatic Filter Media Technology
- Consensus 15: Breathing Resistance



3M Electret Technology Development

- Objective:
 - To develop respirators with high filter performance by advancing the materials and processing of our nonwoven electret technology.
- Benefit:
 - Higher performance filter media will help allow us to reduce respirator thickness, weight, and pressure drop; increase comfort; service life; and expand our design options.



Bioaerosol Interaction with Respirators

- Consensus 1: Safety and Effectiveness
- Consensus 2: Self-Contamination
- Consensus 3: Fomite Transmission
- Consensus 7: Repeated Disinfection Durability



3M Bioaerosol Generating Equipment



Equipment used to evaluate efficacy of antimicrobial treated respirators and compared to *in vitro* test methods

Bioaerosol Test Method Comparison

Conclusions:

- Antimicrobial treatment evaluated on respirator did not appear to reduce the viability of pathogens on the respirator compared to control respirator without antimicrobial treatment
- This was contrary to the outcome with the standard *in vitro* test method

Ylitalo, C., Stepanova, N., Laingen E., Sebastian J., Viner A. Bioaerosol Interaction with Respirators: The Efficacy of Antimicrobial Treatment as tested by a Standard *in vitro* test method and by a bioaerosol test method. Paper presented at American Industrial Hygiene Conference & Expo Session, May 18, 2011



Bioaerosol Filtration

- University of Nebraska Medical Center, Department of Pathology & Microbiology, Bioaerosol Research Laboratory
- 3M Company, Personal Safety Division, St. Paul, MN
- This research was supported by a grant from the Air Force Research Laboratory.

Conclusions:


- Evidence that conventional filtering facepiece respirators can be very effective at capturing virus particles.
- No statistical evidence to demonstrate that an antimicrobial agent on a respirator has an effect on the few virus particles that do penetrate through the filter.

Lore M., Brown T., Hinrichs S., Sebastian J., Viner A., McCullough N. Performance of Conventional & Antimicrobial – Treated Filtering Facepiece Respirators against Viable Influenza A. Paper presented at American Industrial Hygiene Conference & Expo Session PO 124, May 18, 2011.



3M Personal Safety Division

Fluid Resistance - ASTM 1860 Test Apparatus



3M

3M and Project BREATHE

- 3M is a partner with the VA and NIOSH on Project BREATHE
- 3M's contribution to this effort
 - *Fit and Comfort*
 - *Filter Media*
 - *Respirator Design*
 - *Microbial Testing*
 - *Analytical Capability*