

# Filtering Facepiece Respirators and Viable Microbial Aerosols

## NIOSH Personal Protective Technology Program Healthcare Stakeholder Meeting

**Roundtable 3** – Emerging Topic - Considerations for Extending  
Respirator Supplies During an Outbreak or Pandemic

June 18, 2013

**Brian Heimbuch**  
**Senior Bioaerosol Scientist**  
**Associate Division Manager**  
850-914-3188  
bheimbuch@ara.com



# FFR Decontamination Study

N95 Filtering Facepiece Respirators (FFRs) used for Study



3M 1860



3M 1870



3M 8210



3M 8000






Kimberly-Clark




Moldex 2201

# FFR Decontamination Study

Method	Conditions	Graphic
Microwave Generated Steam (MGS)	2-minute cycle on a water reservoir 1250 Watt microwave	
Ultraviolet Germicidal Irradiation (UVGI)	15-minute treatment @ 1.6 – 2.2 mW/cm <sup>2</sup> (1.8 X10 <sup>4</sup> J/M <sup>2</sup> )	
Low-Temperature Moist Heat (WMH)	30 minutes, 65 °C, 85% RH	

# FFR Decontamination Study



**E35 Committee**


- Similar results were obtained using low-path H5N1<sup>1</sup>

**E2721: Standard test method for effectiveness of decontamination of air-permeable materials challenged with biological aerosols containing human pathogenic viruses**

<sup>1</sup>Lore MB, BK Heimbuch, TL Brown, JD Wander, SH Hinrichs, Effectiveness of Three Decontamination Treatments against Influenza Virus Applied to Filtering Facepiece Respirators. *Annals of Occupational Hygiene*, 2012;56(1):92-101



# FFR Decontamination Study



- Similar results were obtained using low-path H5N1<sup>1</sup>

**ASTM 2720: Standard test method for evaluating the effectiveness of decontamination procedures on surfaces challenged with droplets containing human pathogenic viruses**




**E35 Committee**

<sup>1</sup>Lore MB, BK Heimbuch, TL Brown, JD Wander, SH Hinrichs, Effectiveness of Three Decontamination Treatments against Influenza Virus Applied to Filtering Facepiece Respirators. *Annals of Occupational Hygiene*, 2012;56(1):92-101

# FFR Decontamination Study

Arithmetic Means and Standard Deviation MDFF10 Values across the Control and Three Decontamination Conditions for Six FFR Models



\* Indicates a statistically significant reduction ( $P < 0.05$ ) compared with the control

- Fit was not significantly degraded



# FFR Decontamination Study

FFR	UVGI	MGS	WMH	Control
3M 8210	0.41 ± 0.24	0.08 ± 0.03	0.43 ± 0.37	0.62 ± 0.19
3M 8000	1.24 ± 0.22	1.33 ± 0.24	0.70 ± 0.07	0.88 ± 0.12
Moldex 2201	1.26 ± 0.25	1.25 ± 0.31	0.90 ± 0.29	2.05 ± 0.33
KC PFR	1.59 ± 0.27	2.14 ± 0.22	2.16 ± 0.10	2.12 ± 0.41
3M 1870	0.34 ± 0.40	0.52 ± 0.35	1.06 ± 0.56	0.63 ± 0.35
3M 1860S	0.66 ± 0.14	0.98 ± 0.39	0.58 ± 0.07	0.64 ± 0.07

- Filtration efficiency was not significantly degraded



# FFR Cleaning Study

- The Food and Drug Administration (FDA) requires that reprocessed single-use medical devices be cleaned and sterilized, and that their functional performance be demonstrated<sup>2</sup>
- Cleaning studies were performed on 3M1860, 3M1870, and Kimberly-Clark N95 surgical FFRs contaminated with *S. aureus* and artificial saliva using ASTM E2721



Wipe Product	Active Ingredient
Pampers® Wipe	None
3M™ 504/07065 Respirator Cleaning Wipe	Benzalkonium Chloride (BAC)
Current Technology Inc. Hype-Wipe®	Hypochlorite (OCL)

<sup>2</sup>Medical Device User Fee and Modernization Act of 2002, Public Law 107-250.



# FFR Cleaning Study

*S. aureus* Removal using Pampers Wipe    Artificial Saliva Removal using Pampers Wipe



Nose pads


Edge Strips

FFR Fabric

- OCL wipe showed 4 – 5 log reduction of *S. aureus* on most surfaces
  - Nose pads for 3M1870 had 1 – 2 log reduction
- BAC wipe showed 2 – 5 log reduction of *S. aureus* on all surfaces
  - Edge strips of Kimberly-Clark FFR showed 4 – 5 log reduction
- BAC wipe removal of artificial saliva was similar to Pampers wipe

# FFR Cleaning Study

## Effect of Cleaning on Particle Removal Efficiency




# FFR Hospital Wear Study

- Environmental staff at Bay Medical Center (Panama City, FL) wore FFRs while cleaning discharged patient rooms
- FFR wear time was ~25 minutes
- Staff was trained to don and doff the FFRs to avoid contact contamination of the FFRs by touching
- Following doffing, the FFRs were deconstructed and bacterial isolates were collected using permissive media



# FFR Hospital Wear Study

## Viability Bacterial Isolates Recovered From FFRs



- 73% of the Gram-positive and 67% of the Gram-negative isolates evaluated were resistant to oxacillin
- Vancomycin resistance was 9.2% and 36.7%, respectively

# Viable H1N1 Evaluation of FFRs



## Average Removal Efficiencies of 0.8- $\mu$ m Particles at 85 LPM

FFR Model	Inert	H1N1 influenza	<i>p</i> -value
3M 1860S	99.85% $\pm$ 0.10%	99.27% $\pm$ 0.38%	0.08
3M 1870	99.90% $\pm$ 0.09%	99.13% $\pm$ 1.36%	0.45
Kimberly-Clark	99.72% $\pm$ 0.16%	98.93% $\pm$ 0.36%	0.02
SafeLife T5000	99.999% $\pm$ 0.001%	99.996% $\pm$ 0.002% <sup>a</sup>	0.09
GSK Actiprotect	99.94% $\pm$ 0.06%	99.23% $\pm$ 1.00%	0.19

## Average Removal Efficiencies of 0.8- $\mu$ m Particles at 170 LPM

FFR Model	Inert	H1N1 influenza	<i>p</i> -value
3M 1860S	99.37% $\pm$ 0.39%	98.56% $\pm$ 0.87%	0.13
3M 1870	99.96% $\pm$ 0.03%	99.59% $\pm$ 0.27%	0.14
Kimberly-Clark	98.37% $\pm$ 0.32%	96.29% $\pm$ 0.56%	0.02
SafeLife T5000	99.994% $\pm$ 0.009%	99.995% $\pm$ 0.002% <sup>a</sup>	0.9
GSK Actiprotect	99.23% $\pm$ 0.15%	96.29% $\pm$ 2.49%	0.09



# Path Forward

**Data from these studies can be used for both short-term and long-term solutions for mitigating an FFR shortage**

## Short-Term Solution

- FFRs are robust enough to be decontaminated
- FFRs cannot be cleaned according to FDA guidelines, but we do not think cleaning is necessary based on projected operational guidelines
- **Risk cannot be eliminated**, but could be reduced with a little more research aimed at specific risk factors
  1. Strain resistance risk
  2. Repeated exposure complications
  3. Universal application
  4. Increased decontamination cycles
  5. Transition preparations for UVGI technology

# Path Forward

**Data from these studies can be used for both short-term and long-term solutions for mitigating an FFR shortage**

## Long-Term Solution

- Data from these studies can be used to develop better FFRs
  - Reuse is approved within a shift
  - Reuse between patients will be allowed during a pandemic
- Development of a decontaminatable FFR is allowable and could meet current NIOSH and FDA approval guidelines<sup>4</sup>
  - Co-develop cleaning protocols/devices
  - Material and design features can be optimized based on our data

<sup>4</sup>Heimbuch BK, Harnish D. (2011) Discussions on Short-Term and Long-Term Solutions to Mitigate a Shortage of Filtering Facepiece Respirators Caused by Pandemic Influenza, Final Report from Interagency Meeting, Food and Drug Administration-Centers for Devices and Radiologic Health. Available from the Food and Drug Administration

# Acknowledgements

